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
Figure 1. Overview of IPTG induction

1. Prokaryotic cells (*E. coli*) DO NOT produce T7 RNA polymerase.
2. Can genetically insert the following genetic elements (called “DE3”) into prokaryotic chromosome:
 - 甲、 the *lac* operator
 - 乙、 the *lac* promoter
 - 丙、 T7 RNA polymerase gene
3. Under normal circumstance, lac I (produced by the *lac I* gene) will bind to the *lac* operator, thereby preventing the DE3 (which harbors the T7 RNA polymerase gene) from being transcribed by *E. coli* RNA polymerase.
4. Addition of IPTG will replace lac repressor from the DE3, allowing production of T7 RNA polymerase from *E. coli* genome.
5. Now, the newly made T7 RNA polymerase can act on the plasmid carrying the recombinant cDNA (e.g. ECCD23) to induce the production of recombinant protein product.

Figure 2. Maps of pMal-p2x plasmid

Features:

lacI: codes for the *lac* repressor protein

lacZα: codes for galactosidase enzyme

polylinker: equivalent to the “multiple cloning sites”

Amp^r: Ampicillin-resistant gene


Figure 3. A typical IPTG-induced production of ECCD23 fusion proteins



Figure 4. A typical bacterial growth curve

TENS protocol for the Mini-Preps of plasmid DNA

[Procedures]

1. Pellet overnight cultures in microcentrifuge tubes for one minute.
2. Discard supernatant but leave approximately 50 to 100 ul in tube.
3. Add 300 ul of TENS buffer (see below). Vortex.
4. Set on ice for 10 min.
5. Add 150 ul of 3.0 M Sodium Acetate, pH 5.2 and invert.
6. Spin for 10 min.
7. *(Optional) Transfer supernatant to new tube. Extract with 500 µl PCA (Phenol:Chloroform:Iso-amyl alcohol, 25:24:1). Vortex after adding the PCA solution. Spin 13,000 rpm, 10 min.
8. Carefully transfer supernatant to fresh tube containing 0.9 ml of 100% EtOH pre-cooled to -20°C.
9. Spin for 2 min.
10. Discard supernatant. Rinse pellet 1X with 1 ml of 70% EtOH – Spin for 10 min.
11. Discard supernatant. Resuspend pellet with 30 ul of ddH₂O.
12. Determine concentration of the plasmid using spectrophotometer (read OD_{260nm}).

[Materials & Reagents]

3 M Sodium acetate (pH 5.2)	100 ml
10 N NaOH	100 ml
100% EtOH (-20°C)	100 ml
70% EtOH (room temperature)	100 ml
PCA (phenol: chloroform: amyl alcohol=25:24:1)	

TENS buffer

10 mM Tris-HCl, pH 8.0	1 ml of 1 M Tris-HCl, pH 8.0
1 mM EDTA, pH 8.0	200 ul of 0.5 M EDTA, pH 8.0
0.1 N NaOH	1 ml of 10 N NaOH
0.5% SDS	5 ml of 10% SDS
	Add H ₂ O to 100 ml

[Note]

1. (Step #7) The PCA step is optional if the preps are used for restriction digests - however, the plasmid DNA is not likely to be stable over time without PCA extraction due to potential co-purifying DNases.
2. (Step #11) Tris-EDTA (TE) buffer is a better choice for the long-term storage of plasmid DNA at 4C. TE contains EDTA capable of chelating divalent ions (e.g. Mg²⁺) required for the proper function of enzymes (e.g. DNase) detrimental to DNA. But, if the plasmid DNA is meant for the sequencing purpose, then the use of TE buffer MUST BE avoided (use ddH₂O instead).